

## TYGON® tubing LMT-55 - non-DEHP laboratory and vacuum tubing

This tubing is a replacement of the tubing r3603/r3607 by Saint Gobain. The formerly applied DEHP plasticizers are banned in the European market due to EU regulation. The LMT-55 applies allowed plasticizers.

### Features/Benefits

- The new non-DEHP formulation delivers the same excellent performance as the r3603/r3607
- Outstanding chemical resistance
- Lot-to-lot consistency for reproducible results
- Non-oxidizing and non-contaminating
- Smooth, polished inner wall
- Slips easily over fittings and grips securely for simple lab set-ups
- Contains no BPA
- Standard sizes available to hold full vacuum at room temperature

### Specifications

Product type	Tubing
Packaging unit	1 m
Temperature range	-46°C – 74°C
Tubing ID	1.3 mm
Tubing OD	3 mm
Wall thickness	0.85 mm
Sterilizable	Yes (autoclave or ethylene oxide)
Color	Transparent
Manufacturer	Saint Gobain
Certifications	FDA compliant, meets USP Class VI and NSF 51 requirements